	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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Transpro CD 01 plus is a universal culture medium that can be used for the production of adenovirus in suspension of HEK293 cells, as well as for passaging, high-density culture, and transient transfection culture of HEK293 and CHO cells. Transpro CD 01 plus is suitable for the production of adenoviruses, lentiviruses, adeno-associated viruses, etc. during the research and development process using HEK 293 series cells such as HEK 293, Expi293F, 293F, HEK 293-SF-3F6, 293E, etc. Meanwhile, Transpro CD 01 plus is also suitable for transient transfection and expression culture of antibodies and recombinant proteins during the development process using HEK293 series cells such as HEK 293, Expi293F, 293F, 293E, and CHOK1, CHOZN, and other CHO series cells. This product is a fully chemically limited culture medium, free of animal derived ingredients, free of protein components, and free of animal or plant derived protein hydrolysates. Transpro CD 01 plus liquid packaging does not contain L-glutamine, and an additional 4-6 mM L-glutamine needs to be added during use.

【细胞培养 Cell culture】

- ① 建议细胞接种密度 Suggested cell inoculation density:  $0.2\sim1.0\times10^6$  cells/mL.
- ② 温度 Temperature: 37°C
- ③ CO<sub>2</sub>:8%
- ④ 使用时需要额外添加 L-谷氨酰胺 4-6mM

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

It needs to be supplemented with 4-6mM L-glutamine when used.

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Most cells lines can adapt directly into this product. They can be directly inoculated into this medium and passed more than three times. For some cell lines, sequential cell adaptation may be used when using this medium.

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Prepare frozen solution on the super clean workbench : 90% Transpro CD 01 plus +10% dimethyl sulfoxide (DMSO) mixed solution, precooling at 2~8℃ (heat will be released when DMSO is diluted);

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Frozen cell fluid: cells are in the exponential growth period, with a density greater than  $1.5 \times 10^6$  cells/mL, and the viability is greater than 95%.

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Cell fluid was centrifuged at 800rpm for 5 min;

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Slowly pour out the supernatant, resuspend the cells with cryopreservation solution, the cryopreservation density is  $1.0 \sim 1.5 \times 10^7$  cells/mL, and transfer the cells to a sterile cryopreservation tube;

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Place the cryopreservation tube in the cryopreservation box containing isopropyl alcohol, freeze it at - 80 °C overnight, and then transfer it to the liquid nitrogen tank for long-term storage. If there is no freezing box, the temperature can be reduced manually by gradient as follows:

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- -80℃冻存过夜;
- freeze at - 80℃ overnight;
- 转移至液氮罐中长期贮存。
- transfer frozen cells to liquid nitrogen tank for long-term storage.

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Prepare a 37 °C warm water to thaw frozen cells;

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Take out the frozen tube from the liquid nitrogen tank and quickly thaw frozen cells in 37°C warm water;

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After wiping the cryopreservation tube with 75% ethanol, open the cryopreservation tube in the sterile operation table, transfer the cell fluid to a 15 ml centrifuge tube containing 2-5 mL of Transpro CD 01 plus, blow and mix well, centrifuge at 800 rpm for 5 minutes;

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Slowly pour out the supernatant, resuspend with 15~20 ml preheated Transpro CD 01 plus, and transfer to a 125 ml shake flask;

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Place it in a shaking incubator with 8% CO<sub>2</sub>, 80% humidity, 110 ~ 130rpm, at 37°C for culture;

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After 2~3 days of culture, the cells were counted and subcultured.

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The cells are seeded at  $0.2\sim1.0\times10^6$  cells/mL, count and subculture every 2 ~ 3 days. In the first three passages, the volume remained unchanged to restore cell viability. After the cell viability recovers to normal and reaches more than 90%.

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

The seed cells were expanded at the density of  $0.2\sim1.0\times10^6$  cells/mL until reaching the required volume. The criteria for normal seed state: the viability was greater than 95%, the cell morphology was regular and round, and the growth doubling time was normal.

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The day before transfection need to seed cells at  $2.0\times10^6$  viable cells/ml, the cell density can reach  $4.0\times10^6$  viable cells/ml on the second day ;

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After cell counting on the second day of culture, the cell viability was more than 95%, and the living cell density was  $\geq 4.0\times10^6$  cells / ml, can be used directly; If the cell density is lower than  $4.0\times10^6$  cells / ml, the cells can be collected by centrifugation (800 rpm, 5 min), and the cells can be separated at  $4.0\times10^6$  cells / ml density was resuspended in Transpro CD 01 plus medium;

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The mixture of DNA and PEI was prepared according to the optimized transient process;

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Add the mixed solution to the culture medium for culture;

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After 18 hours of culture, it is recommended to supplement the supplemented medium Transpro feed 1 (the concentration is recommended to be 3~5% of the initial culture volume), or the combined supplemented medium DN feed B2 (the concentration is recommended to be 0.3~0.5% of the initial culture volume), which can further improve the density of living cells and protein expression.

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Culture until 7 days, or the vitality is less than 60%, and end the culture.

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Transpro CD 01 plus liquid packaging: 2°C to 8°C, protect from light; validity period: 12 months.

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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