	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

产品使用说明书 Product Use Instructions

多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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Transpro CD 01 (-Gln) medium is a universal transient medium, which can be used for subculture, high-density culture and transient transfection culture of HEK293 cells and CHO cells. The transient transfection process does not require centrifugation to change the medium. Transpro CD 01 (-Gln) medium is suitable for the use of 293 series cells such as HEK 293, Expi293F, 293F, 293E and CHO series cells such as ExpiCHOS and CHOS for transient transfection expression culture of antibodies, recombinant proteins and viruses during the development and manufacture process. Transpro CD 01 (-Gln) medium is an animal-derived component free (ACF), protein free (PF), chemically defined (CD) medium, and does not contain any growth factor and hydrolysates, which ensures consistency between batches and improves the efficiency of the cell culture process. This product does not contain HT, anti-clumping agent or L-Glutamine.

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Suitable for powder package (taking 1L as an example)

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1.Prepare ultrapure water with a volume of about 90% (900 mL);

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2.Add 22.80g of Transpro CD 01 (-Gln) medium powder, and stir for 30min to dissolve completely;

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3.Add sodium bicarbonate 2.220g, and stir for 5~10min until completely dissolved;

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	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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5. Detecting pH and osmotic pressure; (Adjust pH7.00 to 7.40, Osmolality 275 to 320 mOsm/kg)

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6. Sterilize and filter with a 0.22μm filter.

【细胞培养 Cell culture】

① 建议细胞接种密度 suggested cell inoculation density: $0.2\sim1.0\times10^6$ cells/mL

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Most cell lines use this product without any domestication, and can be directly inoculated into this medium and passed for more than three times. For some cell lines, gradient continuous domestication may be used when using this series of medium.

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Prepare the cryopreservation solution on the super clean workbench: 90% Transpro CD 01 + 10% dimethyl sulfoxide (DMSO) mixture, precooling at 2~ 8°C(heat will be released when DMSO is diluted);

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Cryopreserved cell fluid: in logarithmic growth stage, with a density greater than 1.5×10^6 cells / mL, and the activity rate is greater than 95%;

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The cell fluid was centrifuged at 800 rpm for 5 min;

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Slowly pour out the supernatant and resuspend the cells with cryopreservation solution, and the cryopreservation density is $1.0\sim1.5\times10^7$ cells / mL, transfer the cells to the sterile cryopreservation tube;

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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Place the cryopreservation tube in the cryopreservation box containing isopropyl alcohol, freeze it at - 80 °C overnight, and then transfer it to the liquid nitrogen tank for long-term storage. If there is no freezing box, the temperature can be reduced manually by gradient as follows:

- 4℃冻存 30min;
- freeze at 4℃ for 30min;
- -20℃冻存 2~4 小时;
- freeze at -20℃ for 2~4h;
- -80℃冻存过夜;
- freeze at - 80℃ overnight;
- 转移至液氮罐中长期贮存。
- transfer frozen cells to liquid nitrogen tank for long-term storage.

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Prepare a 37 °C warm water to thaw frozen cells ;

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Take out the cryopreservation tube from the liquid nitrogen tank and rapidly thaw frozen cells in a 37℃ water bath;

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After wiping the cryopreservation tube with 75% ethanol, open the cryopreservation tube in the sterile operation table, transfer the cell suspension to a 15 mL centrifuge tube and centrifuge at 800 rpm for 5 min;

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Slowly pour out the supernatant, resuspend it with 15 ~ 20 mL preheated transpro CD 01, and transfer it to a 125 mL shake flask;

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Place it in a shaking incubator with 8% CO₂，110 ~ 130rpm，at 36.5℃for culture;

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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After 2-3 days of culture, the cells are counted and subcultured.

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The cells are seeded at $0.2 \sim 1.0 \times 10^6$ cells/mL, count and subculture every 2 ~ 3 days. In the first three passages, the volume remained unchanged to restore cell viability. After the cell viability recovers to normal and reaches more than 90%. The seed cells were expanded at the density of $0.2 \sim 1.0 \times 10^6$ cells/mL until reaching the required volume. The criteria for normal seed state: the viability was greater than 95%, the cell morphology was regular and round, and the growth doubling time was normal.

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The day before transfection need to seed cells at 2.0×10^6 cells/mL, the cell density can reach 4.0×10^6 cells/mL on the second day;

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After cell counting on the second day of culture, the cell viability was more than 95%, and the living cell density was $\geq 4.0 \times 10^6$ cells / mL, can be used directly; If the cell density is lower than 4.0×10^6 cells / mL, the cells can be collected by centrifugation (800 rpm, 5 min), and the cells can be resuspended at 4.0×10^6 cells/mL density in Transpro CD 01;

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The mixture of DNA and PEI was prepared according to the optimized transient process;

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Add the mixed solution to the medium for culture;

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After 18 hours of culture, it is recommended to supplement the supplemented medium Transpro feed 1 (the concentration is recommended to be 3-5% of the initial culture volume), or the combined supplemented medium DN feed

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

B2 (the concentration is recommended to be 0.3-0.5% of the initial culture volume), which can further improve the density of living cells and protein production.

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Culture until 7 days, or the viability is less than 60%, and end the culture.

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Transpro CD 01 (-Gln) medium Powder packing: 2°C to 8°C, Protect from light; Retest date: 24 months.

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